



**DN-Micro**  
FLUID SOLUTIONS

**CULTURA<sup>®</sup>**



## **CULTURA** - The reliable and convenient way to incubate dipslides for microbial testing

- Supplied with a dipslide tray to hold 18 dipslides
- Temperature range : 25-45°C (Dipslides should be incubated at approx 30°C)
- Supplied with a two year warranty from purchase date
- Weight: 1.1kg
- External Dimension: 310 x 155 x 168
- Internal Dimensions 220 x 120 x 150



**For further information and  
Prices contact DN-Micro**

Unit 315, Tanners Business Centre  
Waterside Mill, Chew Valley Road,  
Greenfield, Oldham, OL3 7NH  
T: 01457 877722  
E: [info@dn-micro.net](mailto:info@dn-micro.net)  
W: [www.dn-micro.net](http://www.dn-micro.net)